



MORE THAN HEAT 30-3000 °C

Material-ID	100037
Hersteller	Dental Direkt
Beschreibung	Reserve
Material Name	DD Nature Zr

Programmnummer	Name	Prozesse				Vortrocknen				Schließen				Vorwärmen				Brennen 1				Brennen 2				Brennen 3				Kühlen			
		Freigeschaltet	Rate [°/min]	Temperatur [°]	Haltezeit [s]	Vakuum	Freigeschaltet	Rate [°/min]	Temperatur [°]	Haltezeit [s]	Vakuum	Freigeschaltet	Rate [°/min]	Temperatur [°]	Haltezeit [s]	Vakuum	Freigeschaltet	Rate [°/min]	Temperatur [°]	Haltezeit [s]	Vakuum	Freigeschaltet	Rate [°/min]	Temperatur [°]	Haltezeit [s]	Vakuum	Freigeschaltet	Rate [°/min]	Temperatur [°]	Haltezeit [s]	Vakuum		
1	Waschbrand					x		440	120							x	45	860	90	x									x	STEP	440	180	
2	Hauptbrand					x		440	240							x	45	830	90	x									x	STEP	440	180	
3	Korrekturbrand					x		440	240							x	45	820	90	x									x	STEP	440	180	
4	Glanzbrand					x		440	120							x	45	815	60										x	STEP	440	180	
5	Glanzbrand für Vollanatomie					x		440	120							x	45	830	60										x	STEP	440	180	

Material-ID	100037
Manufacturer	Dental Direkt
Description	Reserve
Material name	DD Nature Zr

Program number	Name	Processes				Predrying				Closing				Preheating				Firing 1				Firing 2				Firing 3				Cooling			
		Enable	Rate [°/min]	Temperature [°]	Hold time [s]	Vacuum	Enable	Rate [°/min]	Temperature [°]	Hold time [s]	Vacuum	Enable	Rate [°/min]	Temperature [°]	Hold time [s]	Vacuum	Enable	Rate [°/min]	Temperature [°]	Hold time [s]	Vacuum	Enable	Rate [°/min]	Temperature [°]	Hold time [s]	Vacuum	Enable	Rate [°/min]	Temperature [°]	Hold time [s]	Vacuum		
1	Wash firing					x		440	120							x	45	860	90	x									x	STEP	440	180	
2	Main firing					x		440	240							x	45	830	90	x									x	STEP	440	180	
3	Correction firing					x		440	240							x	45	820	90	x									x	STEP	440	180	
4	Glaze firing					x		440	120							x	45	815	60										x	STEP	440	180	
5	Glaze firing for full contours					x		440	120							x	45	830	60										x	STEP	440	180	